

Advanced high throughput capillary plate for protein crystallization

Beat Blattmann, Mareike Kurz and Markus G. Grütter

Institute of Biochemistry, University of Zürich, Winterthurerstrasse 190, CH-8057 Zürich
Paul Reardon, SWISSCI AG Neuhofstrasse 74 CH-6345 Neuheim



Aims:

- capillary crystallization plate based on the counter-diffusion method⁽¹⁾
- applicable for high throughput screening X-tal experiments
- as well as for individually designed X-tal experiments
- compatible with any incubation and imaging system
- easy to work with (loading protein/ solutions, mounting of X-tals).

Solution: Crystalharp

- SBS formatted capillary plate designed for 48 counter diffusion X-tal experiments

- Manufactured and sold by **SWISSCI AG**



A detailed description on plate handling and plate loading can be found on the leaflet (see below).

Polyimide coated quartz capillary


Capillary internal diameter: 0.1 mm
Capillary length on a plate: 200 cm
Capillary length / experiment: 30 mm
V= 240 nl / experiment
Total of <16 µl of protein needed




Crystallization using Crystalharp

Plates tested on different proteins


Crystals were taken to the SLS for diffraction analysis



Capillaries used in a cryo stream - no ice formation




Capillaries are sealed and mounted in a glass sleeve on a standard magnetic CrystalCap.




Cryo solution diffused partly through the capillary, cryo-protected part stays clear, non cryo-protected part turns intransparent on freezing.

Data collection was performed under cryocooling in a N₂-gas stream at the PX(X06SA) beamline at the SLS (Villigen, CH). Using the data collected from crystals grown in capillary the structures of lysozyme⁽²⁾ and AcrB⁽³⁾ could be solved using MR.



AcrB 3.7 Å



lysozyme 1.4 Å

Summary and conclusion

Counter diffusion crystallization of macromolecules in capillary is an easy, cost-effective, and practical procedure for obtaining protein crystals suitable for in situ X-ray data analysis. The counter diffusion process has been used to simultaneously screen for optimal conditions for protein crystal growth, and mix in cryogenic solutions in a single capillary tube. Problems harvesting crystals and difficulties in transportation are reduced to a bare minimum. We can show, that crystals grown in capillary diffract to at least the same resolution as the ones grown by vapour diffusion. A 1.4 Å dataset for lysozyme and a 3.7 Å dataset for AcrB were collected and the structures solved by molecular replacement. Additionally, we observed that capillary grown crystals can be flash-frozen without the need of a cryo protectant. The observation of ice rings was reduced to a bare minimum.

(1) Garcia-Ruiz JM, Methods in Enzymology, 1997, Vol. 368

(2) Cianci M, Helliwell JR, Suzuki A., Acta Crystallogr D, 2008 Dec;64

(3) Pos KM, Schiefner A, Seeger MA, Diederichs K., FEBS Lett. 2004 Apr 30;564(3)